

**Preservative-Free
0.9% Sodium Chloride
Injection, USP**

Rx only

DESCRIPTION

0.9% Sodium Chloride Injection, USP in **Carpject™** sterile cartridge unit (2 mL, 3 mL, and 5 mL fill), is a sterile aqueous isotonic solution of sodium chloride in Water for Injection. The total osmolar concentration of the solution is 0.31 mOsm/mL. Each mL contains 9 mg of sodium chloride.

CLINICAL PHARMACOLOGY

0.9% Sodium Chloride Injection, USP is a sterile aqueous vehicle having approximately the same osmotic pressure and composition as extracellular fluids. It is nonirritating to tissues.

INDICATIONS AND USAGE

For use as a sterile isotonic vehicle for diluting or dissolving compatible parenteral medications.

PRECAUTIONS

If compatibility is in doubt, consult the appropriate specialized literature before mixing with any drug product. Do not administer if a precipitate appears, if solution is cloudy or hazy, or if any other visible change occurs upon mixing.

For single use only. Discard unused solution.

DOSAGE AND ADMINISTRATION

Use in accord with any warnings or precautions appropriate to the medication being administered.

HOW SUPPLIED

NDC No.	Container	Fill	Quantity
0409-1918-32	Carpject™ with Luer Lock	2 mL fill in 2 mL cartridge	bins of 50
0409-1918-33	Carpject™ with Luer Lock	3 mL fill in 5 mL cartridge	bins of 25
0409-1918-35	Carpject™ with Luer Lock	5 mL fill in 5 mL cartridge	bins of 25

Do not use if solution is discolored or contains a precipitate.

Store at 20 to 25°C (68 to 77°F). [See USP Controlled Room Temperature.]

Do not freeze.

Revised: November, 2009

Printed in USA
Hospira, Inc., Lake Forest, IL 60045

EN-2305

